

ABSO-MATE™ PLEATED DEPTH FILTER CARTRIDGE

All polypropylene, absolute-rated, cost-effective filtration

Parker's Abso-Mate™ Cartridges provide the ultimate in economical filtration for even the most critical process fluids. The proprietary melt blown media is rigidly controlled for reliable results time after time. Abso-Mate cartridges are produced without adhesives that can potentially contaminate fluids.


Abso-Mate Pleated Cartridges are available in 0.2µm, 0.45µm, 1µm, 2µm, 5µm, 10µm, 20µm, 40µm, and 70µm absolute rated pore sizes.

Benefits

- Absolute ratings for consistent and reliable performance (99.98%; $\beta = 5000$)
- Back-washable media, reduces replacement maintenance and cartridge disposal costs
- Abso-Mate cartridges are non-fiber releasing and contain minimal extractables
- All materials of construction are FDA listed as acceptable for potable and edible liquid contact according to CFR Title 21
- One-piece construction eliminates bypass concerns on multi-length cartridges
- All-polypropylene construction offers wide chemical compatibility with most chemicals, acids, bases and solvents
- Fused construction and continuous lengths eliminate the need for adhesives and allow accurate bubble point integrity testing
- ISO 9001 registered company

Applications

- Membrane Prefilter
- Chemicals
- Catalyst Recovery
- Precious Metal Recovery
- Waste Water



Specifications

Materials of Construction

Type of Construction: Integrally sealed, all-polypropylene pleated media supported by all-polypropylene construction

Filter Media: Melt blown polypropylene microfiber

Media Support Layers: Non-woven or mesh polypropylene

Media Support Core: Heavy wall high strength polypropylene

Media Support Cage and Thermally Welded End Caps: Molded polypropylene

Seal Materials: Buna-N, EPR, Silicone, Viton®, PFA Encapsulated Viton®

Dimensions

Cartridge Outside Diameter: 2-1¹/₁₆ in.

Cartridge Inside Diameter: DOE: 1-1¹/₁₆ in., SOE: 1-5¹/₃₂ in.

Maximum Recommended Operating Conditions

Temperature: 200°F (93°C)

Change Out ΔP: 35psi (2.4bar)

ΔP @ Ambient 70°F (21°C): 90psi (6bar)

ΔP @ 200°F (93°C): 20psi (1.4bar)

Flow Rate: 10gpm (38 lpm) per 10 in. length

Biological Safety/Product Purity

- All components FDA listed per CFR, Title 21
- Non-fiber releasing per FDA Part 210.3B (5) and (6)
- Non-photo sensitive

Filtration Ratings

99.98% efficiency at 0.2, 0.45, 1, 2, 5, 10, 20, 40, & 70 μm pore sizes.

Beta Ratio (β) =
$$\frac{\text{Upstream Particle Count @ Specified Particle Size and Larger}}{\text{Downstream Particle Count @ Specified Particle Size and Larger}}$$

Percent Removal Efficiency =
$$\left(\frac{\beta-1}{\beta}\right) 100$$

Performance determined per ASTM F-795-88. Single-Pass Test using AC test dust in water at a flow rate of 3.5gpm per 10 in. (13.2 lpm per 254 mm) cartridge.

Performance Attributes

Flow Rate and Pressure Drop Formulas

Flow Rate (gpm) =
$$\frac{\text{Clean } \Delta P \times \text{Length Factor}}{\text{Viscosity} \times \text{Flow Factor}}$$

Clean ΔP =
$$\frac{\text{Flow Rate} \times \text{Viscosity} \times \text{Flow Factor}}{\text{Length Factor}}$$

- Notes:
1. Clean ΔP is psi differential at start.
 2. Viscosity is centistokes. Use Conversion Tables for other units.
 3. Flow Factor is psid/gpm at 1cks for 10 in. (or single).
 4. Length Factors convert flow or ΔP from 10 in. (single length) to required cartridge length.

Abso-Mate Flow Factors
(psid/gpm @ 1 cks)

Rating (μm)	Flow Factor
0.20	3.100
0.45	1.000
1	0.750
2	0.300
5	0.072
10	0.031
20	0.021
40	0.012
70	0.008

Abso-Mate Length Factors

Length (in)	Length Factor
9	1.0
10	1.0
19	2.0
20	2.0
29	3.0
30	3.0
39	4.0
40	4.0

Liquid Particle Retention Ratings (μm) @ Removal Efficiency of:

Cartridge	β=5000 Absolute	β=1000 99.9%	β=100 99%	β=50 98%	β=20 95%
PAB002	0.2	<0.2	<0.2	<0.2	<0.1
PAB004	0.45	0.4	0.2	<0.2	<0.1
PAB010	1	0.8	0.4	<0.2	<0.1
PAB020	2	1.9	0.8	<0.2	<0.1
PAB050	5	3.8	1.4	0.4	0.15
PAB100	10	7	2	0.5	0.25
PAB200	20	13	4	1.8	0.35
PAB400	40	22	7	3.2	0.8
PAB700	70	52	22	15	5.5



ABSO-MATE™ PLEATED DEPTH FILTER CARTRIDGE

Ordering Information



Filter Rating		Nominal Length			Support Construction		Seal Material		End Cap Configuration		Special Options	
CODE	MICRON	CODE	INCHES	MM	CODE	MATERIAL	CODE	MATERIAL	CODE	DESCRIPTION	CODE	DESCRIPTION
002	0.2	9	9 5/8	244	F	Glass-filled Polypropylene (core only)	P	Polyethylene Foam (DOE gasket only)	AR	020 O-ring/Recessed cap	B	Bubble-point test
004	0.45	10	9 13/16	249			E	EPR	DO	Double open end (DOE)	R	DI water rinse (5 min.)
010	1	19	19 5/8	498	G	304 Stainless Steel (core only)	N	Buna-N	DX	Double open end/extended core	Z6	Individual Poly bag only
020	2	20	19 15/16	506			S	Silicone	LL ²	120/120 (Filterlite LMO & Nuclepore Polymeric Vessels)		
050	5	29	29 1/4	743	A	Natural Polypropylene (All support components)	T ¹	PFA Encapsulated Viton® (222, 226 & O-ring only)*	LR ²	120 O-ring/Recessed (Nuclepore)		
100	10	30	30 1/16	764			V	Viton®	OB	Std. open end/Polypropylene spring closed end		
200	20	39	39	991			X	No seal material	PR ²	213 O-ring/Recessed cap (Ametek® & Parker LT Polymeric Vessels)		
400	40	40	40	1016					SC	226 O-ring/Flat		
700	70								SF	226 O-ring/Fin		
									SSC	SS inserted 226 O-ring/Closed		
									SSF	SS inserted 226 O-ring/Fin		
									TC	222 O-ring/Flat		
									TF	222 O-ring/Fin		
									STC	SS inserted 222 O-ring/Closed		
									STF	SS inserted 222 O-ring/Fin		
									TX	222 O-ring/Flex Fin		
									XB	Ext. core open end/Polypropylene spring closed end		

Specifications are subject to change without notification. For User Responsibility Statement, see www.parker.com/safety. Viton is a registered trademark of E.I. DuPont de Nemours & Co., Inc. Ametek is a registered trademark of Ametek, Inc.

Parker Hannifin Corporation
Bioscience & Water Filtration Division
 2340 Eastman Avenue
 Oxnard, California, USA 93030
 +1 805 604 3400
bwf.oxn.support@support.parker.com
www.parker.com/bioscience
www.parker.com/industrial-filtration

DS_ME_Abso-Mate_2025_Rev. A

© 2025 Parker Hannifin Corporation



WARNING: This product can expose you to chemicals including 1,3-Butadiene which are known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

²Available only in 9 5/8" (-9) and 19 5/8" (-10) lengths