

# Connection Solutions for Life Sciences & Clean Room Applications



## Respiratory

Oxygen Therapy/Oxygen Gas/Transfilling/Concentrators/  
Oxygen Conserving Devices/Sleep Apnea/ICU/Aerosol



## Preventative & Recovery

Perioperative Temperature Mgmt./Pre-OP/OR/Post-OP/  
Therapeutic/Compression Therapy/Alternating Support Surfaces



## Bio-Fluid Management

Dialysis/Medical Autoclaves/Dental/Hospital Infection Control/  
Suction Therapy/Wound Therapy

## Surgical & Diagnosis

Surgical Power Tools/Imaging Equipment/Home Diagnostic  
Equipment/Advanced Prosthetics

MEDICAL

INDUSTRIAL

## Laboratory

Gas Control/AA Spectrometry/Thermal Conductivity Detector

## Clean Rooms

Air/Vacuum Conditioning Unit/Air Bearing Controller/  
Semi-Conductor/Neonatal Ventilator/Filling & Packing

## Pharmaceutical

Air & Nitrogen Supply/Buffer Preparation/Bioreactor Production/  
Chromatography/Diafiltration & Concentration/Dosing/  
Filling & Packing

## Respiratory

Anti-Dust Systems/O<sub>2</sub> Delivery Systems

