

CHO-BOND® 1086

Silicone Adhesive Primer

Parker Chomerics CHO-BOND primers 1086 is an air-drying liquid coating used to improve the adhesion of Parker Chomerics CHO-BOND conductive silicone compounds to metal and other non-silicone substrates.

The primer is moisture reactive and clear in color. CHO-BOND 1086 primer is formulated for use with CHO-BOND 1016, 1030, 1035, 1038, 1075 electrically conductive adhesives/sealants and CHO-THERM® 1641 thermal compound.

Note that typically, the proper primer comes bundled as a "kit" with the CHO-BOND or CHO-THERM compound ordered. However, if extra primer is deemed necessary, then part numbers for just the primer are provided in Table 2.

See the Conductive Compounds Selector Guide for a more detailed listing of which CHO-BOND adhesive / sealant part numbers come packaged as a kit with the primer and which do not.

Typical Applications:

For best bonding results, use the following procedure:

1. Surfaces to be joined must be clean, dry and oil free. In a well-ventilated area, clean the substrates with isopropyl alcohol (IPA) and let the solvent flash off. (Other solvents, such as methyl-ethyl-ketone (MEK), toluene, or acetone can be used to clean difficult to remove contaminants).
2. Wet a lint-free cotton cloth with the primer and apply to the surface in horizontal and vertical strokes. Keep the surface wet at all times. A cotton swab is a common applicator for a grooved or stepped surface.

NOTE: Due to the low flash point temperature of the primer, ensure that the container is re-sealed immediately after use.

3. The primed surface must cure for 30 minutes at room temperature. A relative humidity of 40% to 70% is optimal. A low humidity level may require a longer cure time.

NOTE: There is no elevated temperature cure time for the primer.

4. Apply the CHO-BOND/CHO-THERM adhesive/sealant per that product's instructions.

NOTE: If the primer has been applied for more than 4 hours without the application of the adhesive/sealant, repeat steps 1 through 4.



CHO-BOND® 1086 PRODUCT INFORMATION

Typical Properties	CHO-BOND 1086	Test Methods
Viscosity @ 25°C (77°F), centistokes	5	
Color	Clear	
Specific Gravity	0.82	
Flash point, °C (°F)	14 (57)	
Use Temperature, °C (°F)	-80 to 200 (-112 to 392)	
Elevated Temperature Cure	None	
Room Temperature Cure	0.5 hr	
Working Life	N/A	
Shelf Life, unopened, from Date of Manufacture	9 months at 25°C (77°F)	
Minimum thickness recommended	0.0001 in (0.00254 mm)	
Maximum thickness recommended	0.0005 in (0.01270 mm)	
Volatile Organic Content (VOC)	731 g/l	

CHO-BOND® 1086 ORDERING INFORMATION

Product	Weight (grams)	Packaging	Part Number
CHO-BOND 1086	10	3 dram glass vial	50-10-1086-0000
CHO-BOND 1086	95	4 fluid ounce glass bottle	50-04-1086-0000
CHO-BOND 1086	375	1 pint can	50-01-1086-0000



We're Here to Help

Scan QR code or visit parker.com/chomerics to:

- [Request a Free Sample](#)
- [Talk to an Expert](#)
- [Get a Quote](#)
- [Find Where to Buy](#)

Parker Hannifin Corporation

Chomerics Division

77 Dragon Court

Woburn, MA 01801

Phone 781 935 4850

Fax 781 933 4318

chomailbox@parker.com

parker.com/chomerics

CHODS1086 September 2024

©2024 Parker Hannifin Corporation

